

*Safety First  
Healthcare*





Safety First Healthcare developed its business based on in-depth research and established contacts in Asia. When we founded Safety First we decided to investigate Asia countries which could meet four important, basic requirements:

- 1) High government standards and strict regulations for personal protective equipment
- 2) Country free of political entanglements and controversy affecting the supply chain
- 3) Country free of trade wars
- 4) Skilled and reasonably priced labor pool

Based on this criteria, Safety First quickly found a large mask factory in Vietnam, called Nutica, which has Vinacontrol, CE, FDA and ISO certifications, but was primarily selling to Vietnam clients. Safety First is the majority shareholder of Nutica and since taking ownership, we have gathered 22 other small factories in Vietnam under the Nutica umbrella. Together we are able to produce 35 million masks per week and are significantly increasing our production capacity. We have been able to expand the mask product offerings to include a variety of masks specific to the needs of our medical professional clients.

Additionally, Safety First offers masks for the general public and are in the process of expanding the offerings. In this new world environment having comfortable, safe masks suited for all ages has become a critical consideration for most families.

We also have a gown, cap and bootie factory, called Delatech, under the Safety First group. Delatech produces various types of medical gowns as well as full protective and isolation suits with hood and booties. Delatech has Vinacontrol, CE, FDA and ISO certifications.

The pages that follow display our entire suite of personal protective clothing. For further information please contact us at [info@safetyfirst.healthcare](mailto:info@safetyfirst.healthcare).



# MASKS



# NUTICA CERTIFICATE – LICENSE







# PROTECTION FACE MASK

## 4 Ply - Disposable




Safety First  
Healthcare



# PROTECTION FACE MASK

4 Ply - Disposable



● Color:





# PROTECTION FACE MASK

4 Ply – Disposable

## ● Description

This 4 Ply Mask is disposable only. This kind of mask is used for protection of the respiratory system. The mask is hygienic, comfortable and effective. This 4 Ply white antibacterial mask is a breakthrough product of NUT masks. It is more than 95% preventative against small particulates and bacteria. The mask is manufactured from non-woven fabric with a micro-filtering layer and meltblown antibacteria layer that effectively filter at a higher than 95% efficiency. The design has a folding edge around the mask, an elastic band and a nose strap allowing a closer fit. The design makes the mask more airy and easier to breathe and to adapt to many face shapes and dimensions.

## ● Instructions for use

- Do not expose mask to extreme light or heat.
- Do not use if oxygen level is below 19.5%.
- Do not use in environments of harmful substances with high concentration levels.
- Do not use the mask in environments of toxins, amianthus or dust caused by sand blast.

## ● Characteristics

- Elastic band: White elastic material.
- Plastic strap bar (steel covered by PE plastic)
- Outside layer and inside layer (Layer: 1&2): Polypropylen and polyester white or blue non-woven fabric.
- Weight: 30/gsm
- Anti-bacterial layer (Layer 3): Meltblown Paper - SMS. Milky white
- Weight: 25-30/gsm. Thickness: 0.2mm
- Bacteria Filtration Efficiency %: 99%
- Filtration Resistance Pa: 30Pa.
- Outside and Inside layer (Layer 4): Polypropylen and polyester white non-woven fabric. Weight: 20/gms.
- Product color: White or Blue







# PROTECTION FACE MASK

4 Ply - Disposable

## ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

## ● Commitments & Certification of Nutica

## ● Packaging

Piece/box	Box/carton	Piece/carton
50	50	2500

## ● The use events

People who are in any of the below circumstances should wear this mask:

- In crowded environments such as air travel and public transportation.
- When required by law or regulation for contact with individuals who are suspected of having infection, communicable disease, Covid-19 and other contagious illness.
- When contacting or caretaking of individuals with symptoms such as cough, runny nose, difficulty breathing or other respiratory tract illnesses.

Safety First  
Healthcare



**NUT**<sup>®</sup>  
Nutica

# PROTECTION FACE MASK

## 3 PLY DISPOSABLE



## ● Description


This 3 Ply Mask is disposable only. This kind of mask is used for protection of the respiratory system. The mask is hygienic, comfortable and effective. This 3 Ply white antibacterial mask is a breakthrough product of NUT masks. It is more than 95% preventative against small particulates and bacteria. The mask is manufactured from non-woven fabric with a micro-filtering layer and meltblown antibacteria layer that effectively filter at a higher than 95% efficiency. The design has a folding edge around the mask, an elastic band and a nose strap allowing a closer fit. The design makes the mask more airy and easier to breathe and to adapt to many face shapes and dimensions.

## ● Instructions for use

- Do not expose mask to extreme light or heat.
- Do not use if oxygen level is below 19.5%.
- Do not use in environments of harmful substances with high concentration levels.
- Do not use the mask in environments of toxins, amianthus or dust caused by sand blast.

## ● Characteristics

- Elastic band: White elastic material.
- Plastic strap bar (steel covered by PE plastic)
- Outside layer (Layer 1): Polypropylen and polyester white or blue non-woven fabric.
- Weight: 30/gsm
- Anti-bacterial layer (Layer 2): Meltblown Paper - SMS. Milky white
- Weight: 25-30/gsm. Thickness: 0.2mm
- Bacteria Filtration Efficiency %: 99%
- Filtration Resistance Pa: 30Pa.
- Outside and Inside layer (Layer 3): Polypropylen and polyester white non-woven fabric. Weight: 20/gms.
- Product color: White or Blue





# PROTECTION FACE MASK

## 3 PLY DISPOSABLE

### ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

### ● Commitments & Certification of Nutica

### ● Packaging

Piece/box	Box/carton	Piece/carton
50	50	2500

### ● The use events

People who are in any of the below circumstances should wear this mask:

- In crowded environments such as air travel and public transportation.
- When required by law or regulation for contact with individuals who are suspected of having infection, communicable disease, Covid-19 and other contagious illness.
- When contacting or caretaking of individuals with symptoms such as cough, runny nose, difficulty breathing or other respiratory tract illnesses.



# 3D MEDICAL MASK

SUPER PROTECTION



Outstanding advantages:  
Breathable  
Comfortable - ready fit  
Does not wrinkle by skin  
Does not cause discomfort behind ears due to long straps  
Design to limit fogging of the glasses when wearing




100% CERTIFIED  
NUT-10000  
EXP-10000



# 3D MEDICAL MASK

SUPER PROTECTION



## ● Description

Mask is disposable only. Mask is used to protect the respiratory system in an effective, comfortable and hygienic manner. 4-layer white disposable mask for protection from dust and droplets. It is manufactured from a non-woven fabric including the micro-filtering layer. The product is a breakthrough of NUT masks. It prevents small particle size penetration and is antibacterial. The mask achieves a filtering efficiency higher than 95%. The mask has a 3D folding edge and a elastic wearing band. 3D mask is designed in a parallel piped shape to give the nose room and to adapt to the many shapes and dimensions of various faces. It also helps prevent fogging of glasses.

## ● Instructions for use

- Do not expose mask to extreme light or heat.
- Do not use if oxygen level is below 19.5%.
- Do not use in environments of harmful substances with high concentration levels.
- Do not use the mask in environments of toxins, amianthus or dust caused by sand blast.

## ● Characteristics

- Color: White.
- Elastic band: White elasticity substance.
- Non-woven fabric and milky white antibacterial fabric.
  - Outside layer (Layer 1): Polypropylen and polyester non-woven fabric. Weight: 30/gms.
  - Inside layer (Layer 2): Polypropylen and polyester non-woven fabric. Weight: 25/gms.
  - Antibacterial layer (Layer 3): Spunbond + Meltbown - SMS + Spunbond Nonwoven fabric. Weight: 25-30/gms.
- Inside layer (Layer 4): Polypropylen and polyester non-woven fabric. Weight: 25/gms.





# 3D MEDICAL MASK

SUPER PROTECTION

## ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

## ● Commitments & Certification of Nutica

## ● Packaging

Piece/box	Box/carton	Piece/carton
10	50	500

## ● The use events

Recommendation: The people who are in the below circumstances should wear the mask:

- When coming to crowded places or travel on by public means.
- When obligatorily required for contact, looking after of a person who is infected or a person who is suspect of being infected from Virus Corona epidemic disease.
- Should wear a mask when looking after or contacting people who have symptoms of respiratory tract disease, such as, having a cough, difficulty breathing, runny nose and so on.
- For isolation or disease examination or disease treatment at disease examination or treatment establishments.




Safety First  
Healthcare



NUT<sup>®</sup>  
Nutica

# MEDICAL FACE MASK NUT95





## ● Description

NUT 95 Respirators. NUT 95 respirators are respiratory protective devices designed to achieve a very snug fit and efficient filtration of airborne particles. The 'NUT 95' designation means that when subjected to careful testing, the respirator will block at least 95 percent of very small (0.3 micron) test particles.

## ● Instructions for use

NUT 95 respirators are used as protection by medical personnel in the operating room during surgical procedures. The respirators protect the surgical operator and the surgical patient from the transfer of body fluids, particulate pollution, and micro organisms.

## ● Characteristics

- Layer 1: Polypropylene (SS) Non-woven (50g)
  - Layer 2: Polypropylene (PP) 95 melt-blown nonwovens(25g)
  - Layer 3: Polypropylene (PP) 95 melt-blown nonwovens(25g)
  - Layer 4: Tuyere cotton(50g)
  - Layer 5: Polypropylene (SS) Non-woven (20~30g)
- At the point of nose folding is made by polyethylene, steel band inside.  
Chinlon is the elastic material.





# NUT95

## ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

## ● Commitments & Certification of Nutica

## ● Packaging

Piece/box	Box/carton	Piece/carton
28	50	1400

## ● The use events

People who are in any of the below circumstances should wear this mask:

- In crowded environments such as air travel and public transportation.
- When required by law or regulation for contact with individuals who are suspected of having infection, communicable disease, Covid-19 and other contagious illness.
- When contacting or caretaking of individuals with symptoms such as cough, runny nose, difficulty breathing or other respiratory tract illnesses.

Safety First  
Healthcare



NUT<sup>®</sup>  
Nutica

# 4PLY MASK KIDS



Safety First  
Healthcare



NUT<sup>®</sup>  
Nutica

# 4PLY MASK KIDS



Safety First  
Healthcare



# 4PLY MASK KIDS





# 4PLY MASK KIDS

## Instruction for use

- Do not expose mask to extreme light or heat.
- Do not use if oxygen level is below 19.5%.
- Do not use in environments of harmful substances with high concentration levels.
- Do not use the mask in environments of toxins, amianthus or dust caused by sand blast.

## Warranty

The sole responsibility of the seller or manufacturer shall replace the quantity of product which has been demonstrated to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as, the returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not and the user must bear every responsibility and risk relating to such use.

## Certificate, license

This product is certified by the National Institute of Medical Device and Construction under the Ministry of Health of Vietnam, No. 082820/VTTB-DGCL in accordance with the technical specifications of Standard TCVN 8389-1:2010- Normal medical mask. To be licensed for use to protect the respiratory system, dust filtering, antibacterial and water-proofing. Pursuant to Basic Standard of Nutica Company, No. TCCS: 01-2020/CTY-NUTICA and certificate, FDA – CE – ISO 9001:2015 – ISO 13485:2016 (Further supplementing the issued licenses).

## The use events

People who are in any of the below circumstances should wear this mask:

- In crowded environments such as air travel and public transportation.
- When required by law or regulation for contact with individuals who are suspected of having infection, communicable disease, Covid-19 and other contagious illness.
- When contacting or caretaking of individuals with symptoms such as cough, runny nose, difficulty breathing or other respiratory tract illnesses.

## Commitments & Certification of Nutica







# 4PLY MASK KIDS

## Description

This 4-Ply Children's Mask is disposable only. This kind of mask is used for protection of the respiratory system. The mask is hygienic, comfortable and effective for children. This 4-ply blue and white antibacterial mask is a breakthrough product of NUT masks. It is 95% preventative against small particulates and bacteria. The mask is manufactured from non-woven fabric with a micro-filtering layer that effectively filters at a higher than 95% efficiency. The design has a folding edge around the mask, an elastic band and a nose strap allowing a closer fit. The design makes the mask more airy and easier to breathe and to adapt to many face shapes and dimensions.

## Characteristics

- Elastic band: White elastic material.
- Plastic strap bar (steel covered by PE plastic)
- Outside layer and inside layer  
(Layer: 1&2): Polypropylen and polyester white or blue non-woven fabric.
- Weight: 30/gsm
- Anti-bacterial layer (Layer 3):  
Spunbond + Meltblown + Spunbond Nonwoven Fabric Milky
- Weight: 25-30/gsm.
- Outside and Inside layer (Layer 4): Polypropylen and polyester white or blue non-woven fabric.
- Weight: 20/gms.
- Product is blue or white

## Product composition:

- Elastic string form Polyester àn spandex
- Nonwoven fabric Polypropylene Spunbond SS
- Nonwoven fabric Polypropulene Meltblown
- Plastic HD nasal splint bar




## Packaging

Piece / Box	Box / Carton	Piece / Carton
50	50	2500

Safety First  
Healthcare



# 3D MASK KIDS



Safety First  
Healthcare




# 3D MASK KIDS





# 3D MASK KIDS





# 3D MASK KIDS





# 3D MASK KIDS



## ● Description

Mask is disposable only. Mask is used to protect the respiratory system in an effective, comfortable and hygienic manner. 4-layer white disposable mask is for protection from dust and droplets. It is manufactured from a non-woven fabric including the micro-filtering layer. The product is a breakthrough of NUT masks. It prevents small particle size penetration and is antibacterial. The mask achieves a filtering efficiency higher than 95%. The mask has a 3D folding edge and an elastic ear strap band. 3D mask is designed in a parallel piped shape to give the nose room and to adapt to the many shapes and dimensions of various faces. It also helps prevent fogging of glasses.

## ● Instructions for use

- Do not expose mask to extreme light or heat.
- Do not use if oxygen level is below 19.5%.
- Do not use in environments of harmful substances with high concentration levels.
- Do not use the mask in environments of toxins, amianthus or dust caused by sand blast.

## ● Characteristics

### Composition:

- Soft wearing band, non-woven SMS Plus 3+ and SS waterproof

### Usage and advantages:

- Preventing dust, bacteria and the communicable diseases from the respiratory tract
- Prevents UV rays from skin
- 3D design for a close face fit and airy design making it easier to breathe





# 3D MASK KIDS

## ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

## ● Commitments & Certification of Nutica

## ● Packaging

Piece/box	Box/carton	Piece/carton
25	100	2500

## ● The use events

People who are in any of the below circumstances should wear this mask:

- In crowded environments such as air travel and public transportation.
- When required by law or regulation for contact with individuals who are suspected of having infection, communicable disease, Covid-19 and other contagious illness.
- When contacting or caretaking of individuals with symptoms such as cough, runny nose, difficulty breathing or other respiratory tract illnesses.



# **PROTECTIVE CLOTHING**



# DELATECH CERTIFICATE - LICENSE





**Vinacontrol**  
Since 1977

**CERTIFICATE OF ANALYSIS**

No. 000170/2016/01

Item	REQUIREMENT	Result
1	Identification of the product	OK
2	Traceability of the product	OK
3	Control of the product	OK
4	Storage of the product	OK
5	Transportation of the product	OK
6	Delivery of the product	OK

**Remarks:**

1. The number of samples is consistent.  
2. The test equipment is calibrated.  
3. The test method is correct.  
4. The test results are consistent with the requirements.  
5. The test results are consistent with the requirements.

**Technical Manager:** VINACONTROL - TESTING LAB (VINACONTROL S1)  
**Tester:** TRAN THI HOANG  
**Reviewer:** TRAN THI HOANG

Page 1/1

**Vinacontrol**  
Since 1977

**CERTIFICATE OF ANALYSIS**

No. 000170/2016/01

Item	REQUIREMENT	Result
1	Identification of the product	OK
2	Traceability of the product	OK
3	Control of the product	OK
4	Storage of the product	OK
5	Transportation of the product	OK
6	Delivery of the product	OK

**Remarks:**

1. The number of samples is consistent.  
2. The test equipment is calibrated.  
3. The test method is correct.  
4. The test results are consistent with the requirements.  
5. The test results are consistent with the requirements.

**Technical Manager:** VINACONTROL - TESTING LAB (VINACONTROL S1)  
**Tester:** TRAN THI HOANG  
**Reviewer:** TRAN THI HOANG

Page 1/1

**Vinacontrol**  
Since 1977

**PHIẾU KẾT QUẢ THỬ NGHIỆM**  
Certificate of analysis

No. 000170/2016/01

Item	REQUIREMENT	Result
1	Identification of the product	OK
2	Traceability of the product	OK
3	Control of the product	OK
4	Storage of the product	OK
5	Transportation of the product	OK
6	Delivery of the product	OK

**Remarks:**

1. The number of samples is consistent.  
2. The test equipment is calibrated.  
3. The test method is correct.  
4. The test results are consistent with the requirements.  
5. The test results are consistent with the requirements.

**Technical Manager:** VINACONTROL - TESTING LAB (VINACONTROL S1)  
**Tester:** TRAN THI HOANG  
**Reviewer:** TRAN THI HOANG

Page 1/1

**Vinacontrol**  
Since 1977

**PHIẾU KẾT QUẢ THỬ NGHIỆM**  
Certificate of analysis

No. 000170/2016/01

Item	REQUIREMENT	Result
1	Identification of the product	OK
2	Traceability of the product	OK
3	Control of the product	OK
4	Storage of the product	OK
5	Transportation of the product	OK
6	Delivery of the product	OK

**Remarks:**

1. The number of samples is consistent.  
2. The test equipment is calibrated.  
3. The test method is correct.  
4. The test results are consistent with the requirements.  
5. The test results are consistent with the requirements.

**Technical Manager:** VINACONTROL - TESTING LAB (VINACONTROL S1)  
**Tester:** TRAN THI HOANG  
**Reviewer:** TRAN THI HOANG

Page 1/1



# PERSONAL PROTECTION KITS



# PERSONAL PROTECTION KITS

## ● Description

Personal Protective Kit consists of the following:

- Coveralls with cap
- Medical face mask
- Medical Glasses/face shield
- Medical gloves
- Shoe covers

Materials are supplied by large-scale manufacturing company. The manufacturing procedures require is that quality is checked from the input materials to the output products. These products are environmentally friendly and do not causes negative effects to the user's health or environment. These products are onetime use, disposable only.

## ● Instructions for use

- For the coveralls you lightly place your feet into legs then insert your hands into sleeves and put the cap on your head and pull up the zipper. The next will be shoe covers, mask, gloves, glasses/face shield
- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purposes.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.





# PERSONAL PROTECTION KITS

## ● Main materials

- Coveralls with cap: 100% Nonwoven
- Medical face mask: : 100% Nonwoven
- Medical Glasses/face shield: 100% PET
- Medical gloves: 100% rubber(Nitrile rubber/Natural rubber latex)
- Shoe covers: 100% Nonwoven

## ● Warranty

The sole responsibility of the seller or manufacturer shall be to replace defective product if caused by the factory. In the event of product defects, the customer may call our customer service line at +84 906 869 772. Customer service shall notify the customer of the settlement process depending on each specific event (such as, the returned product, refund money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or a consequence of using this product. Before use, the use is requested to determine whether the product is appropriate for the purpose it will be worn or not and the user must bear all responsibility and risk relating to such use.

## ● The use events

- For wearing clean protective clothing when entering a patient's room or at the patient care area
- This is disposable protective clothing and must be eliminated after use.



# MEDICAL PROTECTIVE CLOTHING



# NON - WOVEN FABRIC PROTECTIVE SUIT

## ● Description

Protective clothing is sewed with non-woven fabric, using an elastic and thread combination. This protective clothing is manufactured from the Polypropylene (PP) non-woven fabric. Materials are supplied by large-scale woven fabric manufacturing companies. The manufacturing procedures require that quality is checked from the input materials to the output product. This product is environmentally friendly and does not cause negative effects to the user's health or environment. The product is onetime use, disposable only.

## ● Instructions for use

- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purpose.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.

## ● Main materials

- 100% Non - wovens fabric







# NON - WOVEN FABRIC PROTECTIVE SUIT

## ● Warranty

The sole responsibility of the seller or manufacturer shall replace the quantity of product which has been demonstrated to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement depending on each specific event (such as the returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not and the user must bear every responsibility and risk relating to such use.

## ● Commitments & Certification

## ● Vacuum Packaging

GSM	Qty/bag	Qty/carton	G.W (approx)
30	1	100	13.8
40	1	90	16.0
50	1	80	17.4
60	1	70	18.0
Level 4 suit 60 GSM $\pm$ 5	1	60	16.0
Full Cover body suit 70 GSM	1	40	16.0

## ● The use events

- For wearing clean protective clothing when entering a patient's room or at the patient care area
- This is disposable protective clothing and must be eliminated after use.

Safety First  
Healthcare



# MEDICAL PROTECTIVE GOWN



# NON - WOVEN FABRIC PROTECTIVE GOWN

## ● Description

Protective gown is sewed with non-woven fabric, using an elastic and thread combination. This protective clothing is manufactured from the Polypropylene (PP) non-woven fabric. Materials are supplied by large-scale woven fabric manufacturing companies. The manufacturing procedures require that quality is checked from the input materials to the output product. This product is environmentally friendly and does not cause negative effects to the user's health or environment. The product is onetime use, disposable only.

## ● Instructions for use

- Lightly insert arms into sleeves to wear the coat, tie the band at the neck first, then tie the band in the back.
- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purpose.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.

## ● Main materials

- 100% Non - wovens fabric





# NON - WOVEN FABRIC PROTECTIVE GOWN

## ● Warranty

The sole responsibility of the seller or manufacturer will be to replace the quantity of products which have been found to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer will not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not, and the user must bear all responsibility and risk relating to such use.

## ● Commitments & Certification

## ● Vacuum Packaging

GSM	Qty/bag	Qty/carton	G.W (approx)
20	1	150	12.2
25	1	150	13.7
30	1	150	15.2
40	1	130	16.7
Level 4 gown 60 GSM $\pm$ 5	1	100	16.6

## ● The use events

- For wearing clean protective clothing when entering a patient's room or at the patient care area
- This is disposable protective clothing and must be eliminated after use.



# MEDICAL CAP



## ● Description

Protective cap is sewed with non-woven fabric, using an elastic and thread combination. This protective cap is manufactured from the Polypropylene (PP) non-woven fabric. Materials are supplied by large-scale woven fabric manufacturing companies. The manufacturing procedures require that quality is checked from the input materials to the output product. This product is environmentally friendly and does not cause negative effects to the user's health or environment. The product is onetime use, disposable only.

## ● Instructions for use

- Lightly put over the head and make sure hair is tucked in
- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purpose.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.

## ● The use events

- For wearing clean protective clothing when entering a patient's room or at the patient care area
- This is disposable protective clothing and must be eliminated after use.

## ● Warranty

The sole responsibility of the seller or manufacturer shall replace the quantity of product which has been demonstrated to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement method depending on each specific event (such as the returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not and the user must bear every responsibility and risk relating to such use.

## ● Main materials

- 100% Non-woven fabric


## ● Packaging

- 2000 Pcs/carton





# MEDICAL SHOE COVER



# MEDICAL SHOE COVER

## ● Description

Protective socks is sewed with non-woven fabric, using an elastic and thread combination. This protective socks is manufactured from the Polypropylene (PP) non-woven fabric. Materials are supplied by large-scale woven fabric manufacturing companies. The manufacturing procedures require that quality is checked from the input materials to the output product. This product is environmentally friendly and does not cause negative effects to the user's health or environment. The product is onetime use, disposable only.

## ● Main materials

- 100% Non - wovens fabric

## ● Packaging

- 2000 Pcs/Carton

## ● Warranty

The sole responsibility of the seller or manufacturer shall replace the quantity of product which has been demonstrated to be defective at the factory. Please contact customer service at this number (+84 28 3636 8837) with any concerns. Customer service shall notify the customer of the settlement mehto depending on each specific event (such as the returned product, refunded money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or as a consequence of using this product. Before use, the user is requested to determine whether the product is appropriate with the purpose for its use or not and the user must bear every responsibility and risk relating to such use.

## ● Instructions for use

- Lightly wrap around feet
- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purpose.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.

## ● The use events

- For wearing clean protective clothing when entering a patient's room or at the patient care area
- This is disposable protective clothing and must be eliminated after use.







# PROTECTIVE GLASSES



## ● Description

Modern eye protection bringing together safety, functionality, quality and design, keeps the eyes safe from particulates, water or chemicals

- Technical Parameters
- Material : Polycarbonate, PP
- Lens color : Clear
- Frame color : Black
- Frame size : Resizable

Gender : Unisex

## ● Instructions for use

- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purpose.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.

## ● Main materials

- 100% PET





# PROTECTIVE GLASSES

## ● Warranty

The sole responsibility of the seller or manufacturer shall be to replace defective product if caused by the factory. In the event of product defects, the customer may call our customer service line at (+84 2839633448). Customer service shall notify the customer of the settlement process depending on each specific event (such as, the returned product, refund money, replacement for another product and so on). Both the seller and manufacture shall not be responsible for any personal injury, loss or damage, whether it is direct or an consequence of using this product. Before use, the use is requested to determine whether the product is appropriate for the purpose it will be worn or not and the user must bear all responsibility and risk relating to such use.

## ● Packaging

Qty/carton	Carton/size	G.W (approx)
200	29/43/65cm	7.7kg

## ● The use events

People who are in any of the below circumstances should wear this mask:

- Cycling, Driving, Golf, Hiking, Mountain Biking, Running, Tactical, Triathlon.



# FACE SHIELD



## ● Description

MEDICAL PROTECTIVE FACE SHIELD is a part of Personal Protective Kit. that are used by many users for protection of the facial area and associated mucous membranes (eyes, nose, mouth) from splashes, sprays, and spatter of body fluids. This product is made with PET and thread combination. Materials are supplied by large-scale manufacturing company. The manufacturing procedures require is that quality is checked from the input materials to the output products. These products are environmentally friendly and do not causes negative effects to the user's health or environment. These products are onetime use, disposable only.

## ● Main Materials

100% PET

## ● Technical Parameter,

- Size: 30cm Width x 22cm height or customerized
- Gender: Unisex

## ● Warranty

The sole responsibility of the seller or manufacturer shall be to replace defective product if caused by the factory. In the event of product defects, the customer may call our customer service line at +84 906 869 772. Customer service shall notify the customer of the settlement process depending on each specific event (such as, the returned product, refund money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or a consequence of using this product. Before use, the use is requested to determine whether the product is appropriate for the purpose it will be worn or not and the user must bear all responsibility and risk relating to such use.

## ● Packaging

Qty/bag	Qty/carton	Carton size	G. W (approx)
01	500	62/35/75 cm	19 kg



# FACE SHIELD




## ● Description

MEDICAL PROTECTIVE FACE SHIELD is a part of Personal Protective Kit. that are used by many users for protection of the facial area and associated mucous membranes (eyes, nose, mouth) from splashes, sprays, and spatter of body fluids. This product is made with PET and thread combination. Materials are supplied by large-scale manufacturing company. The manufacturing procedures require is that quality is checked from the input materials to the output products. These products are environmentally friendly and do not causes negative effects to the user's health or environment. These products are onetime use, disposable only.

## ● Main Materials

100% PET





# FACE SHIELD

## ● Warranty

The sole responsibility of the seller or manufacturer shall be to replace defective product if caused by the factory. In the event of product defects, the customer may call our customer service line at +84 906 869 772. Customer service shall notify the customer of the settlement process depending on each specific event (such as, the returned product, refund money, replacement for another product and so on). Both the seller and manufacturer shall not be responsible for any personal injury, loss or damage, whether it is direct or a consequence of using this product. Before use, the use is requested to determine whether the product is appropriate for the purpose it will be worn or not and the user must bear all responsibility and risk relating to such use.

## ● Technical Parameter,

- Size: 30cm Width x 22cm height or customized
- Gender: Unisex
- Non-sterile
- Onetime use only

## ● Packaging

Qty/bag	Qty/carton	Carton size	G. W (approx)
01	500	62/35/75 cm	19 kg

## ● Instructions For Use

- Wear the face shield on your head to cover the entirety of user's face.
- This product is used for personal protection, water-proofing, protection from dust and for anti-bacterial purposes.
- Keep product in a dry, clean place, keep away from water and direct contact with the chemical substances to avoid damaging product quality.